

· Bandages specialists since 1956 ·

Calvo Izquierdo designs and manufactures bandages and medical textiles. Founded in 1956, we are a family owned company located in Valencia, Spain.

Our goal is to market high quality products that fulfill the needs of healthcare professionals and final users alike. In order to assure continued product excellence, we have assembled an outstanding team and employ modern facilities and production technologies.

Our portfolio includes products such as cohesive bandages, crepe bandages, paper bandages and elastic bandages.

In addition to technological products related to the health field.

Specially designed device for measuring bandage compression

Our goal is to provide the exact pressure measurement of your compression bandages. With this system, you will be able to measure the pressure of your bandages at a specific moment or over a concrete time period. The pressure is measured by 4 sensors and the results are shown in touch screen as decimal numbers or as a graphic. In both case like static and dynamic conditions. The information is saved in a SD card, storing time and date of the results, with the option of transferring the data to a computer. The Primere + has two different programs depending on the objective: Teaching program and Research Program.



C/ De la Forja 4, (Pol. Ind. El Romeral) Requena 46340 (Valencia) SPAIN

Primere +

What pressure do your bandages have?



Main features

Four Sensors

Four input sensor with more than one meter of length, to measure four points simultaneously.

Store data

Long term data record. Include SD card for store the information with the time and the date.

Two programs

Primere + has two different programs depending on the objectives: Teaching program and research program.

Multiplatform

Easy use touchscreen with differ ent languages.





What pressure do your bandages have?

Developed in collaboration with healthcare students and professionals of the sector

The Primere + has two different programs depending on the objective:

Teaching program: This software has been developed in

·Mode 1: Shows the measurement periodically for each time period selected (between 1 and 9 seconds).

dynamic one.

- ·Mode 2: Shows the measurement every time the "Measure" button is pressed.
- ·Mode 3: Shows graphically the evolution of the more if a proper use is done. measurements.

Data results can be saved and stored in a SD card. This option is available for the three modes mentioned before.

Research program: This software has been developed in order to help students while they are applying different order to help health professionals to monitor bandages' compression bandages techniques. It has three modes in compression meanwhile the practical application. It order to measure the pressure, in a static situation or in a measures the pressure of the bandage periodically in intervals of 1 minute, 15 minutes or one hour. Then, it shows the measurement on the screen and stores the values on the SD card, so the evolution of the bandage can be checked at the end process.

> The useful life of the Primere's + battery is around 36 hours according to the Research program. But it can last